

TiBase, Abutment Screw AND ScanPost

ORDER FORM

Laboratory

Customer ID

Name

First name

Street, number

City, postal code

Phone, Fax

E-mail



IMPLANT MANUFACTURER/SYSTEM			TiBase ¹⁾		Abutment Screw ²⁾		ScanPost ³⁾		
Implant diameter	Platform		Quantity	REF	Quantity	REF	Quantity	REF	Size
Biomet 3i (External hex connection)	3.4	3.4	TiBase B 0 3.4 L	6282557	Abutment Screw B 0 3.4, 4.1, 5.0	6460468	ScanPost B 0 3.4 L	6431089	L
	4.1	4.1	TiBase B 0 4.1 L	6282565			ScanPost B 0 4.1 L	6431105	L
	5	5.0	TiBase B 0 5.0 L	6282573			ScanPost B 0 5.0 L	6431113	L
Biomet 3i Certain® (Internal Connection)	3.4	3.4	TiBase B C 3.4 S	6308048	Abutment Screw B C 3.4, 4.1, 5.0	6460450	ScanPost B C 3.4 S	6431212	S
	4.1	4.1	TiBase B C 4.1 L	6308097			ScanPost B C 4.1 L	6431220	L
	5	5.0	TiBase B C 5.0 L	6308121			ScanPost B C 5.0 L	6431238	L
Dentsply Implants Astra Tech OsseoSpeed TX	3.5 S / 4.0 S	3.5 / 4.0	TiBase AT OS 3.5/4.0 L	6282532	Abutment Screw AT OS 3.5/4.0	6460344	ScanPost AT OS 3.5/4.0 L	6431055	L
	4.5 / 5.0 / 5.0 S	4.5 / 5.0	TiBase AT OS 4.5/5.0 L	6282540	Abutment Screw AT OS 4.5/5.0	6460443	ScanPost AT OS 4.5/5.0 L	6431063	L
Dentsply Implants Frialit / Xive	3.4	3.4	TiBase FX 3.4 S	6282433	Abutment Screw FX 3.4, 3.8, 4.5, 5.5	6460476	ScanPost FX 3.4 S	6430891	S
	3.8	3.8	TiBase FX 3.8 S	6282441			ScanPost FX 3.8 S	6430909	S
	4.5	4.5	TiBase FX 4.5 L	6282458			ScanPost FX 4.5 L	6430917	L
	5.5	5.5	TiBase FX 5.5 L	6282466			ScanPost FX 5.5 L	6430925	L

IMPLANT MANUFACTURER/SYSTEM			TiBase ¹⁾		Abutment Screw ²⁾		ScanPost ³⁾		Size
Implant diameter	Platform		Quantity	REF	Quantity	REF	Quantity	REF	
Nobel Biocare Replace (tri-channel internal connection)	3.5	NP	TiBase NB RS 3.5 L	6282474	Abutment Screw NB RS 3.5	6460526	ScanPost NB RS 3.5 L	6430933	L
	4.3	RP	TiBase NB RS 4.3 L	6282482	Abutment Screw NB RS 4.3, 5.0, 6.0	6460534	ScanPost NB RS 4.3 L	6430941	L
	5	WP	TiBase NB RS 5.0 L	6282490			ScanPost NB RS 5.0 L	6430958	L
	6	6.0	TiBase NB RS 6.0 L	6282508			ScanPost NB RS 6.0 L	6430982	L
Nobel Biocare Nobel Active (internal conical connection)	3.5	NP	TiBase NB A 4.5 L	6308188	Abutment Screw NB A 4.5	6460484	ScanPost NB A 4.5 L	6431279	L
	4.3 / 5.0	RP	TiBase NB A 5.0 L	6308253	Abutment Screw NB A 5.0	6460492	ScanPost NB A 5.0 L	6431287	L
Nobel Biocare Branemark® (External hex connection)	3.3	NP	TiBase NB B 3.4 L	6282516	Abutment Screw NB B 3.4	6460500	ScanPost NB B 3.4 L	6431006	L
	3.75 / 4.0	RP	TiBase NB B 4.1 L	6282524	Abutment Screw NB B 4.1	6460518	ScanPost NB B 4.1 L	6431022	L
Straumann Standard (Tissue Level)	3.3	NN (3.5 mm)	TiBase SSO 3.5 L	6284231	Abutment Screw SSO 3.5	6460559	ScanPost SSO 3.5 L	6431162	L
	3.3 / 4.1 / 4.8	RN (4.8 mm)	TiBase SSO 4.8 L	6284249	Abutment Screw SSO 4.8, 6.5	6460567	ScanPost SSO 4.8 L	6431170	L
	4.8	WN (6.5 mm)	TiBase SSO 6.5 L	6284256			ScanPost SSO 6.5 L	6431196	L
Straumann Bone Level	3.3	NC (3.3 mm)	TiBase S BL 3.3 L	6308154	Abutment Screw S BL 3.3, 4.1	6460542	ScanPost S BL 3.3 L	6431246	L
	4.1 / 4.8	RC (4.1 mm / 4.8 mm)	TiBase S BL 4.1 L	6308337			ScanPost S BL 4.1 L	6431253	L
Zimmer Tapered Screw-Vent	3.7 / 4.1	3.5	TiBase Z TSV 3.5 L	6282581	Abutment Screw Z TSV 3.5, 4.5, 5.7	6460575	ScanPost Z TSV 3.5 L	6431139	L
	4.7	4.5	TiBase Z TSV 4.5 L	6282599			ScanPost Z TSV 4.5 L	6431147	L
	6	5.7	TiBase Z TSV 5.7 L	6282607			ScanPost Z TSV 5.7 L	6431154	L
Medentika M-Implant	3.5 / 4.0 / 4.5 / 5.0	3.5/5.0	TiBase MI 3.5 / 5.0 L	6308295	not available				L

THE FOLLOWING COMPONENTS ARE AVAILABLE, DEPENDING ON THE CONNECTION:

SIZE	Scanbodies for Omnicam		Scanbodies for Bluecam		inCoris ZI meso F0.5		inCoris ZI meso F2	
	REF	Quantity	REF	Quantity	REF	Quantity	REF	Quantity
S	6431311		6431295		6231802		6231828	
L	6431329		6431303		6231810		6231836	

¹⁾ 1x titanium base, 1x Abutment Screw, 1x Scanbody for Bluecam, ²⁾ 2x Abutment Screw, ³⁾ 1x ScanPost, 1x Abutment Screw